

For Research Use Only

Grancalcin Polyclonal antibody

Catalog Number: 15787-1-AP



Basic Information

Catalog Number:

15787-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8486

GenBank Accession Number:

BC005214

GeneID (NCBI):

25801

UNIPROT ID:

P28676

Full Name:

grancalcin, EF-hand calcium binding protein

Calculated MW:

217 aa, 24 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : MCF-7 cells, HL-60 cells, human placenta tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15787-1-AP (Grancalcin antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.