

For Research Use Only

# Grancalcin Polyclonal antibody

Catalog Number: 15787-1-AP



## Basic Information

Catalog Number:

15787-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8486

GenBank Accession Number:

BC005214

GeneID (NCBI):

25801

UNIPROT ID:

P28676

Full Name:

grancalcin, EF-hand calcium binding protein

Calculated MW:

217 aa, 24 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : MCF-7 cells, HL-60 cells, human placenta tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15787-1-AP (Grancalcin antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.