

For Research Use Only

# UBE1L Polyclonal antibody

Catalog Number: 15818-1-AP

2 Publications



## Basic Information

Catalog Number:

15818-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8692

GenBank Accession Number:

BC006378

GeneID (NCBI):

7318

UNIPROT ID:

P41226

Full Name:

ubiquitin-like modifier activating enzyme 7

Calculated MW:

1012 aa, 112 kDa

Observed MW:

112 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse liver tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jie Jin	34257278	Cell Death Dis	WB
Mei-Juan Zhou	28035359	Int J Mol Med	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 15818-1-AP (UBE1L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.