

For Research Use Only

RNASE1 Polyclonal antibody

Catalog Number: 15892-1-AP

1 Publications



Basic Information

Catalog Number:

15892-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8679

GenBank Accession Number:

BC005324

GeneID (NCBI):

6035

UNIPROT ID:

P07998

Full Name:

ribonuclease, RNase A family, 1 (pancreatic)

Calculated MW:

156 aa, 18 kDa

Observed MW:

18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Mouse, rat, human

Cited Species:

mouse

Positive Controls:

WB : mouse pancreas tissue, mouse stomach tissue, rat pancreas tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yong Jin	39205658	Am J Physiol Renal Physiol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15892-1-AP (RNASE1 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.