

For Research Use Only

# CALML4 Polyclonal antibody

Catalog Number: 15894-1-AP

Featured Product

3 Publications



## Basic Information

Catalog Number:

15894-1-AP

Size:

1600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8705

GenBank Accession Number:

BC009516

GeneID (NCBI):

91860

UNIPROT ID:

Q96GE6

Full Name:

calmodulin-like 4

Calculated MW:

196 aa, 22 kDa

Observed MW:

22-25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

## Applications

Tested Applications:

IP, WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, rat

Positive Controls:

WB : mouse kidney tissue, mouse thymus tissue, human colon tissue, rat kidney tissue

IP : mouse thymus tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Myoung Soo Choi	32209652	J Biol Chem	WB,IF
Lian Liu	34448569	ACS Chem Neurosci	WB
Alexandra Holló	37690683	J Biol Chem	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15894-1-AP (CALML4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



IP result of anti-CALML4 (IP:15894-1-AP, 4ug; Detection:15894-1-AP 1:300) with mouse thymus tissue lysate 2800ug.