

For Research Use Only

# CALML4 Polyclonal antibody

Catalog Number: 15894-1-AP

Featured Product

3 Publications



## Basic Information

### Catalog Number:

15894-1-AP

### Size:

1600 µg/ml

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG8705

### GenBank Accession Number:

BC009516

### GeneID (NCBI):

91860

### UNIPROT ID:

Q96GE6

### Full Name:

calmodulin-like 4

### Calculated MW:

196 aa, 22 kDa

### Observed MW:

22-25 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB 1:2000-1:10000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

## Applications

### Tested Applications:

IP, WB, ELISA

### Cited Applications:

WB, IF

### Species Specificity:

human, mouse, rat

### Cited Species:

human, rat

### Positive Controls:

WB : mouse kidney tissue, mouse thymus tissue, human colon tissue, rat kidney tissue

IP : mouse thymus tissue,

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Myoung Soo Choi	32209652	J Biol Chem	WB, IF
Lian Liu	34448569	ACS Chem Neurosci	WB
Alexandra Holló	37690683	J Biol Chem	WB

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15894-1-AP (CALML4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



IP result of anti-CALML4 (IP:15894-1-AP, 4ug; Detection:15894-1-AP 1:300) with mouse thymus tissue lysate 2800ug.