

For Research Use Only

# NECAP2 Polyclonal antibody

Catalog Number: 15899-1-AP

1 Publications



## Basic Information

Catalog Number:

15899-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8721

GenBank Accession Number:

BC018914

GeneID (NCBI):

55707

UNIPROT ID:

Q9NVZ3

Full Name:

NECAP endocytosis associated 2

Calculated MW:

263 aa, 28 kDa

Observed MW:

35-37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : MCF7 cells, HeLa cells, human liver tissue, mouse cerebellum tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhichao Lu	36531205	Oxid Med Cell Longev	WB,IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



MCF7 cells were subjected to SDS PAGE followed by western blot with 15899-1-AP (NECAP2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.