

For Research Use Only

NECAP2 Polyclonal antibody

Catalog Number: 15899-1-AP **1 Publications**



Basic Information

Catalog Number:

15899-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8721

GenBank Accession Number:

BC018914

GeneID (NCBI):

55707

UNIPROT ID:

Q9NVZ3

Full Name:

NECAP endocytosis associated 2

Calculated MW:

263 aa, 28 kDa

Observed MW:

35-37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : MCF7 cells, HeLa cells, human liver tissue, mouse cerebellum tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhichao Lu	36531205	Oxid Med Cell Longev	WB,IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF7 cells were subjected to SDS PAGE followed by western blot with 15899-1-AP (NECAP2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.