

For Research Use Only

QRSL1 Polyclonal antibody

Catalog Number: 15926-1-AP



Basic Information

Catalog Number:

15926-1-AP

Size:

220 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8756

GenBank Accession Number:

BC014389

GeneID (NCBI):

55278

UNIPROT ID:

Q9HOR6

Full Name:

glutaminyl-tRNA synthase
(glutamine-hydrolyzing)-like 1

Calculated MW:

528 aa, 57 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: human placenta tissue, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 15926-1-AP (QRSL1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.