

For Research Use Only

UBR4 Polyclonal antibody

Catalog Number: 15966-1-AP



Basic Information

Catalog Number:

15966-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8740

GenBank Accession Number:

BC018694

GeneID (NCBI):

23352

UNIPROT ID:

Q5T4S7

Full Name:

ubiquitin protein ligase E3

component n-recognin 4

Calculated MW:

5183 aa, 570 kDa

Observed MW:

570 kDa, 270 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IP 0.5-4.0 µg for 1.0-3.0 mg of total

protein lysate

IHC 1:200-1:800

IF 1:50-1:500

Applications

Tested Applications:

IF/ICC, IHC, IP, WB, ELISA

Species Specificity:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: HEK-293 cells, HeLa cells, mouse kidney tissue

IP: HeLa cells,

IHC: human pancreas cancer tissue, mouse testis tissue, human lung cancer tissue

IF: HeLa cells,

Background Information

UBR4 (ubiquitin N-recognin domain-containing E3 ligase 4), also named P600, is an E3 ligase that can bind N-end rule substrates. UBR4-deficient mice die during embryogenesis and exhibit pleiotropic abnormalities (PMID: 23431188). UBR4 plays a role in membrane morphogenesis and extracellular matrix-mediated cell survival and signaling, human papillomavirus (HPV) E7-mediated cellular immortalization, and migration and development of neurons (PMID: 33658508, 25084275). The full-length UBR4 protein has a predicted size of ~570 kDa.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 15966-1-AP (UBR4 antibody) at dilution of 1:400 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 15966-1-AP (UBR4 antibody) at dilution of 1:400 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Various lysates were subjected to Tris-acetate gel system followed by western blot with 15966-1-AP (UBR4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using UBR4 antibody (15966-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



IP result of anti-UBR4 (IP:15966-1-AP, 4ug; Detection:15966-1-AP 1:2000) with HeLa cells lysate 1760 ug.