

For Research Use Only

PAPSS2 Polyclonal antibody

Catalog Number: 15978-1-AP



Basic Information

Catalog Number:

15978-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8766

GenBank Accession Number:

BC009894

GeneID (NCBI):

9060

UNIPROT ID:

O95340

Full Name:

3'-phosphoadenosine 5'-phosphosulfate synthase 2

Calculated MW:

614 aa, 70 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HepG2 cells, human liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 15978-1-AP (PAPSS2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.