

For Research Use Only

SOCS2 Polyclonal antibody

Catalog Number: 16038-1-AP



Basic Information

Catalog Number: 16038-1-AP	GenBank Accession Number: BC010399	Purification Method: Antigen affinity purification
Size: 240 µg/ml	GeneID (NCBI): 8835	
Source: Rabbit	UNIPROT ID: O14508	
Isotype: IgG	Full Name: suppressor of cytokine signaling 2	
Immunogen Catalog Number: AG8909	Calculated MW: 198 aa, 22 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xinyue Gu	38145323	FASEB J	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

