

For Research Use Only

BCCIP Polyclonal antibody

Catalog Number: 16043-1-AP

Featured Product

7 Publications



Basic Information

Catalog Number:

16043-1-AP

Size:

600 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8934

GenBank Accession Number:

BC009771

GeneID (NCBI):

56647

UNIPROT ID:

Q9P287

Full Name:

BRCA2 and CDKN1A interacting

protein

Calculated MW:

314 aa, 36 kDa

Observed MW:

50-58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 ug for 1.0-3.0 mg of total
protein lysate

IF/ICC 1:400-1:1600

Applications

Tested Applications:

WB, IF/ICC, IP, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells, A549 cells

IP : HeLa cells,

IF/ICC : HeLa cells,

Background Information

BCCIP (or TOK1) is isolated as a BRCA2- and CDKN1A (p21)-interacting protein and is implicated in homologous recombination (HR) pathway and inhibition of DNA replication stress. BCCIP is an evolutionarily conserved nuclear protein with multiple interacting domains. BCCIP deficiency in mice impaired embryonic and postnatal neural development, causing severe ataxia, cerebral and cerebellar defects, and microcephaly, which are associated with spontaneous DNA damage and subsequent cell death in the proliferative cell populations of the neural system during embryogenesis. BCCIP is essential for maintaining the transactivation activity of wild type p53 suggesting a potential role of BCCIP in cancer etiology.

Notable Publications

Author	Pubmed ID	Journal	Application
Yi Sui	29932276	FEBS J	WB,IF
Fang Wang	32736693	Biochem Biophys Res Commun	WB
Jiaming Su	27535137	Protein Cell	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 16043-1-AP (BCCIP antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



IP result of anti-BCCIP (IP:16043-1-AP, 4ug; Detection:16043-1-AP 1:800) with HeLa cells lysate 3200ug.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using BCCIP antibody (16043-1-AP) at dilution of 1:800 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).