

For Research Use Only

SAP30L Polyclonal antibody

Catalog Number: 16056-1-AP



Basic Information

Catalog Number: 16056-1-AP	GenBank Accession Number: BC009829	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 79685	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q9HAJ7	
Isotype: IgG	Full Name: SAP30-like	
Immunogen Catalog Number: AG9024	Calculated MW: 183 aa, 21 kDa	
	Observed MW: 21-25 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse, rat

Positive Controls:

WB : human peripheral blood platelets, mouse lung tissue, rat lung tissue

Background Information

SAP30L, also named as NS4ATP2, is one of the key proteins in a multi - subunit protein complex involved in transcriptional regulation via histone deacetylation. SAP30L is an essential component of the Sin3A corepressor complex. The MW migrates to 26-30 kDa with phospho and SUMO modification and may generate complexes.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16056-1-AP (SAP30L antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.