

For Research Use Only

TrpRS Polyclonal antibody, PBS Only

Catalog Number: 16081-1-PBS

Featured Product



Basic Information

Catalog Number:

16081-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9079

GenBank Accession Number:

BC017489

GeneID (NCBI):

7453

UNIPROT ID:

P23381

Full Name:

tryptophanyl-tRNA synthetase

Calculated MW:

471 aa, 53 kDa

Observed MW:

53 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, FC (Intra), IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

TrpRS is also named as IFI53, WRS, WARS and belongs to the class-I aminoacyl-tRNA synthetase family. It catalyzes activation of tryptophan by ATP and transfer to tRNA^{Trp}, ensuring translation of the genetic code for tryptophan. TrpRS is thought to be more closely related to TyrRS than to any other class I enzyme, but the significance of this relationship has been questioned (PMID:7743129). It can exist as a dimer (PMID:7743129). It has 2 isoforms produced by alternative splicing.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 16081-1-AP (TrpRS antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.




Various lysates were subjected to SDS PAGE followed by western blot with 16081-1-AP (TrpRS antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 16081-1-AP (TrpRS antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 16081-1-AP (TrpRS antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.




Various lysates were subjected to SDS PAGE followed by western blot with 16081-1-AP (TrpRS antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.




Various lysates were subjected to SDS PAGE followed by western blot with 16081-1-AP (TrpRS antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.



IP result of anti-TrpRS (IP:16081-1-AP, 4ug; Detection:16081-1-AP 1:500) with HepG2 cells lysate 3000ug. This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed HepG2 cells using TrpRS antibody (16081-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.



1×10^6 HepG2 cells were intracellularly stained with 0.25 ug TrpRS Polyclonal antibody (16081-1-AP) and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Rabbit IgG control Rabbit PolyAb (30000-0-AP) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 16081-1-PBS in a different storage buffer formulation.