

For Research Use Only

TCTN3 Polyclonal antibody

Catalog Number: 16085-1-AP **2 Publications**



Basic Information

Catalog Number: 16085-1-AP	GenBank Accession Number: BC009494	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 26123	Recommended Dilutions: WB 1:200-1:1000 IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q6NUS6	
Isotype: IgG	Full Name: tectonic family member 3	
Immunogen Catalog Number: AG9084	Calculated MW: 607 aa, 66 kDa	
	Observed MW: 66 kDa	

Applications

Tested Applications: IHC, WB, ELISA	Positive Controls: WB : mouse liver tissue, IHC : human lung tissue,
Cited Applications: IF, WB	
Species Specificity: human, mouse	
Cited Species: human, mouse	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

TCTN3 (tectonic-3 or TECT3) is a type I membrane protein that belongs to the tectonic family, which consists of TCTN1, TCTN2 and TCTN3. Studies in mice suggest that tectonic may be involved in Hedgehog (Hh) signaling, and essential for ciliogenesis. Tectonic is component of the tectonic-like complex, a complex localized at the transition zone of primary cilia and acting as a barrier that prevents diffusion of transmembrane proteins between the cilia and plasma membranes. TCTN3 is necessary for transduction of the sonic hedgehog (SHH) signaling pathway. Orofaciodigital syndromes (OFDSs) consist of a group of heterogeneous disorders characterized by abnormalities in the oral cavity, face, and digits and associated phenotypic abnormalities. Recently, TCTN3 mutations are found as the cause of an extreme form of OFD associated with bone dysplasia, tibial defect, cystic kidneys, and brain anomalies.

Notable Publications

Author	Pubmed ID	Journal	Application
Antonia Wiegering	29650680	EMBO J	IF
Garcia-Gonzalo Francesc R FR	21725307	Nat Genet	WB,IF


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung using 16085-1-AP (TCTN3 antibody) at dilution of 1:50 (under 40x lens).



mouse liver tissue were subjected to SDS PAGE followed by western blot with 16085-1-AP (TCTN3 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.