

For Research Use Only

P4HTM Polyclonal antibody

Catalog Number: 16089-1-AP



Basic Information

Catalog Number: 16089-1-AP	GenBank Accession Number: BC011710	Purification Method: Antigen affinity purification
Size: 150 µg/ml	GeneID (NCBI): 54681	
Source: Rabbit	UNIPROT ID: Q9NXG6	
Isotype: IgG	Full Name: prolyl 4-hydroxylase, transmembrane (endoplasmic reticulum)	
Immunogen Catalog Number: AG9123	Calculated MW: 563 aa, 63 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data