

For Research Use Only

GNA11 Polyclonal antibody

Catalog Number: 16103-1-AP



Basic Information

Catalog Number: 16103-1-AP	GenBank Accession Number: BC001528	Purification Method: Antigen affinity purification
Size: 190 µg/ml	GeneID (NCBI): 2767	
Source: Rabbit	UNIPROT ID: P29992	
Isotype: IgG	Full Name: guanine nucleotide binding protein (G protein), alpha 11 (Gq class)	
Immunogen Catalog Number: AG9002	Calculated MW: 170 aa, 19 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data