

For Research Use Only

PHYHD1 Polyclonal antibody

Catalog Number: 16125-1-AP



Basic Information

Catalog Number:

16125-1-AP

Size:

267 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9108

GenBank Accession Number:

BC009373

GeneID (NCBI):

254295

UNIPROT ID:

Q5SRE7

Full Name:

phytanoyl-CoA dioxygenase domain
containing 1

Calculated MW:

297 aa, 33 kDa

Observed MW:

30-32 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human liver tissue, LO2 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 16125-1-AP (PHYHD1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.