

For Research Use Only

ZMIZ2-Specific Polyclonal antibody

Catalog Number: 16158-1-AP



Basic Information

Catalog Number:

16158-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_031449

GeneID (NCBI):

83637

UNIPROT ID:

Q8NF64

Full Name:

zinc finger, MIZ-type containing 2

Calculated MW:

97 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

ZMIZ2 and ZMIZ1 members of a PIAS -like family of proteins that interact with nuclear hormone receptors. This antibody is specific for ZMIZ2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

