

For Research Use Only

GAGE-family Polyclonal antibody

Catalog Number: 16170-1-AP



Basic Information

Catalog Number:

16170-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_021123

GeneID (NCBI):

2579

UNIPROT ID:

POCL80

Full Name:

G antigen 7

Calculated MW:

13 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

GAGE cancer/testis antigens are frequently expressed in many different types of cancer, whereas their expression in normal tissues is limited to the germ cells of the immune-privileged organs, testis and ovary. The GAGE proteins are very close in sequence. This antibody can bind GAGE 1 to -12.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

