

GMEB1-Specific Polyclonal antibody

Catalog Number: 16189-1-AP

Basic Information

Catalog Number:

16189-1-AP

Size:

30 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006582

GeneID (NCBI):

10691

UNIPROT ID:

Q9Y692

Full Name:

glucocorticoid modulatory element
binding protein 1

Calculated MW:

63 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

GMEB1 is a trans-acting factor that binds to glucocorticoid modulatory elements (GME), acting in estrogen receptor-regulated pathway GREB1 is expressed in proliferating prostate tissue & prostate cancer and is thought to function in hormone-dependent cancer growth. It associates with GMEB2 protein, and the complex is essential for parvovirus DNA replication. The protein acts complicatedly in SDS-PAGE and may have multiple electrophoresis forms.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

