

For Research Use Only

# Collagen Type XVIII Polyclonal antibody



Catalog Number: 16191-1-AP

1 Publications

## Basic Information

Catalog Number:

16191-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_030582

GeneID (NCBI):

80781

UNIPROT ID:

P39060

Full Name:

collagen, type XVIII, alpha 1

Calculated MW:

178 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IF, IHC

Species Specificity:

human

Cited Species:

human, rat

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Hao Zhou	38515609	Bioact Mater	IHC, IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

