

PLEKHB2 Polyclonal antibody

Catalog Number: 16265-1-AP

Basic Information

Catalog Number:

16265-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9330

GenBank Accession Number:

BC013991

GeneID (NCBI):

55041

UNIPROT ID:

Q96CS7

Full Name:

pleckstrin homology domain
containing, family B (evectins)
member 2

Calculated MW:

222 aa, 25 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

PLEKHB2, also named as EVT2 (Evectin-2), is involved in retrograde transport of recycling endosomes. PLEKHB2 has been identified as post-Golgi proteins of unknown function that have a PH domain that typically binds polyphosphoinositides. PLEKHB2 is expressed in a broad range of tissues. PLEKHB2 has some isoforms with MW 19-26 kDa and 34 kDa. This antibody recognizes all the isoforms of PLEKHB2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

