

For Research Use Only

OCC1 Polyclonal antibody

Catalog Number: 16306-1-AP



Basic Information

Catalog Number:

16306-1-AP

Size:

170 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9286

GenBank Accession Number:

BC013920

GeneID (NCBI):

387882

UNIPROT ID:

Q8TAD7

Full Name:

chromosome 12 open reading frame

75

Calculated MW:

80 aa, 9 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Background Information

OCC1, also named as C12orf75 and AGD3, is a gene whose expression is particularly abundant in neurons in the macaque primary visual cortex (V1).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

