

For Research Use Only

MGAT4B Polyclonal antibody

Catalog Number: 16308-1-AP



Basic Information

Catalog Number:

16308-1-AP

Size:

70 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9293

GenBank Accession Number:

BC009464

GeneID (NCBI):

11282

UNIPROT ID:

Q9UQ53

Full Name:

mannosyl (alpha-1,3-)-glycoprotein

beta-1,4-N-

acetylglucosaminyltransferase,

isozyme B

Calculated MW:

548 aa, 63 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

