

For Research Use Only

GFOD2 Polyclonal antibody

Catalog Number: 16340-1-AP



Basic Information

Catalog Number:

16340-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9518

GenBank Accession Number:

BC000757

GeneID (NCBI):

81577

UNIPROT ID:

Q3B7J2

Full Name:

glucose-fructose oxidoreductase
domain containing 2

Calculated MW:

42 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 16340-1-AP (GFOD2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.