

For Research Use Only

MRPS21 Polyclonal antibody, PBS Only

Catalog Number: 16359-1-PBS



Basic Information

Catalog Number: 16359-1-PBS	GenBank Accession Number: BC004566	Purification Method: Antigen affinity purification
Concentration: 1 mg/ml	GeneID (NCBI): 54460	
Source: Rabbit	UNIPROT ID: P82921	
Isotype: IgG	Full Name: mitochondrial ribosomal protein S21	
Immunogen Catalog Number: AG8343	Calculated MW: 87 aa, 11 kDa	

Applications

Tested Applications:
IHC, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas tissue slide using 16359-1-AP (MRPS21 antibody) at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 16359-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using MRPS21 antibody (16359-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 16359-1-PBS in a different storage buffer formulation.