

For Research Use Only

ALX1 Polyclonal antibody

Catalog Number: 16372-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

16372-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9295

GenBank Accession Number:

BC010923

GeneID (NCBI):

8092

UNIPROT ID:

Q15699

Full Name:

ALX homeobox 1

Calculated MW:

326 aa, 37 kDa

Observed MW:

35-37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse embryo tissue, ROS1728 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Paul P R Iyyanar	35127681	Front Cell Dev Biol	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse embryo tissue were subjected to SDS PAGE followed by western blot with 16372-1-AP (ALX1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.