

For Research Use Only

NFX1 Polyclonal antibody

Catalog Number:16387-1-AP



Basic Information

Catalog Number:

16387-1-AP

Concentration:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9448

GenBank Accession Number:

BC012151

GeneID (NCBI):

4799

UNIPROT ID:

Q12986

Full Name:

nuclear transcription factor, X-box binding 1

Calculated MW:

1120 aa, 124 kDa

Observed MW:

124 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, Jurkat cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16387-1-AP (NFX1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.