

For Research Use Only

PI4KAP2 Polyclonal antibody

Catalog Number: 16422-1-AP **1 Publications**



Basic Information

Catalog Number:

16422-1-AP

Size:

90 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9619

GenBank Accession Number:

BC020225

GeneID (NCBI):

375133

Full Name:

phosphatidylinositol 4-kinase,
catalytic, alpha pseudogene 2

Calculated MW:

539aa,61 kDa; 577aa,66 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Safiyah Ziyad	29386109	Cell Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data