

For Research Use Only

Caveolin-1 Polyclonal antibody

Catalog Number: 16446-1-AP **3 Publications**



Basic Information

Catalog Number:

16446-1-AP

Concentration:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10924

GenBank Accession Number:

BC006432

GeneID (NCBI):

857

UNIPROT ID:

Q03135

Full Name:

caveolin 1, caveolae protein, 22kDa

Calculated MW:

22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiansheng Ye	32167773	J Agric Food Chem	WB
Jorge L Jimenez-Macias	40009687	Sci Adv	WB
Xiaohui Sun	37716034	Phytomedicine	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data