

For Research Use Only

INTS12 Polyclonal antibody, PBS Only

Catalog Number:16455-1-PBS



Basic Information

Catalog Number:

16455-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9403

GenBank Accession Number:

BC014442

GeneID (NCBI):

57117

UNIPROT ID:

Q96CB8

Full Name:

integrator complex subunit 12

Calculated MW:

462 aa, 49 kDa

Observed MW:

49 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

INTS12 is a subunit of the Integrator complex that is thought to act as a reader to link the complex to chromatin.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 16455-1-AP (INTS12 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16455-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 16455-1-AP (INTS12 antibody) at dilution of 1:500 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 16455-1-PBS in a different storage buffer formulation.