

For Research Use Only

# PCSK5 Polyclonal antibody

Catalog Number: 16470-1-AP

Featured Product

4 Publications



## Basic Information

Catalog Number:

16470-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9552

GenBank Accession Number:

BC012064

GeneID (NCBI):

5125

UNIPROT ID:

Q92824

Full Name:

proprotein convertase  
subtilisin/kexin type 5

Calculated MW:

102 kDa, 207 kDa

Observed MW:

210 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

## Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IHC, IF

Species Specificity:

human, mouse

Cited Species:

human, rat

**Note-IHC: suggested antigen retrieval with  
TE buffer pH 9.0; (\*) Alternatively, antigen  
retrieval may be performed with citrate  
buffer pH 6.0**

Positive Controls:

IHC : human kidney tissue, human brain tissue, human  
skin tissue, human heart tissue, human testis tissue,  
human spleen tissue, mouse brain tissue

## Background Information

PC5, a mammalian kex2-like processing endoprotease, mediates posttranslational endoproteolytic processing for several integrin alpha subunits. It has been identified in rat and mouse endocrine and neuroendocrine tissues, and is encoded by multiple mRNAs that generate two splice variant isoforms referred to as PC5A (102KDa) and PC5B (207KDa). PC5A is smaller than PC5B and it can be processed into two forms (117 kDa and 65 kDa) in neuronal cells (PMID:11457520).

## Notable Publications

Author	Pubmed ID	Journal	Application
Yuki Tachida	28955916	Biochem Biophys Rep	
Jie Tong	33313955	Mol Med Rep	
Mengda Xu	38910562	Circ Heart Fail	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 16470-1-AP (PCSK5 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 16470-1-AP (PCSK5 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).