

For Research Use Only

CD302 Polyclonal antibody

Catalog Number: 16472-1-AP



Basic Information

Catalog Number:

16472-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9557

GenBank Accession Number:

BC020646

GeneID (NCBI):

9936

UNIPROT ID:

Q8IX05

Full Name:

CD302 molecule

Calculated MW:

170aa, 19 kDa; 232aa, 26 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HL-60 cells, U937 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16472-1-AP (CD302 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.