

For Research Use Only

# SNAPC2 Polyclonal antibody

Catalog Number: 16506-1-AP



## Basic Information

Catalog Number:

16506-1-AP

Size:

170 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9689

GenBank Accession Number:

BC011868

GeneID (NCBI):

6618

UNIPROT ID:

Q13487

Full Name:

small nuclear RNA activating complex, polypeptide 2, 45kDa

Calculated MW:

334 aa, 36 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

## Background Information

SNAPC2, also named as SNAP45, is a subunit of SNAPc complex. It recruits TBP and BRF2.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

