

For Research Use Only

VHL Polyclonal antibody

Catalog Number: 16538-1-AP **5 Publications**



Basic Information

Catalog Number: 16538-1-AP	GenBank Accession Number: BC058831	Purification Method: Antigen affinity purification
Size: 650 µg/ml	GeneID (NCBI): 7428	Recommended Dilutions: IHC 1:50-1:500 IF 1:10-1:100
Source: Rabbit	UNIPROT ID: P40337	
Isotype: IgG	Full Name: von Hippel-Lindau tumor suppressor	
Immunogen Catalog Number: AG9843	Calculated MW: 172 aa, 20 kDa	

Applications

Tested Applications: IF/ICC, IHC, ELISA	Positive Controls:
Cited Applications: IF, IHC	IHC: human liver cancer tissue, human kidney tissue
Species Specificity: human	IF: HepG2 cells,
Cited Species: human, mouse	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

VHL (von Hippel-Lindau tumor suppressor) gene was identified as the tumor suppressor gene whose germ line mutations are associated with the inherited von Hippel-Lindau cancer syndrome (PMID: 8603073, 10722748). VHL patients develop a wide variety of tumors including retinal angioma, central nervous system hemangioblastoma, pheochromocytoma, and renal clear cell carcinoma (PMID: 10722748). VHL localizes predominantly to the cytoplasmic compartment but engages in a dynamic nuclear-cytoplasmic shuttle (PMID: 8700833, 9891082, 12101228). VHL has 3 isoforms with the molecular mass of 24, 20 and 18 kDa.

Notable Publications

Author	Pubmed ID	Journal	Application
Fabao Liu	26309161	Oncotarget	
Rui Zhu	34760693	Front Oncol	IHC
Changhui Yu	35351157	Respir Res	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 16538-1-AP (VHL antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of HepG2 cells, using VHL antibody 16538-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).