

For Research Use Only

GSTA4-specific Polyclonal antibody

Catalog Number: 16547-1-AP



Basic Information

Catalog Number:
16547-1-AP

Size:
N/A

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
NM_001512

GeneID (NCBI):
2941

UNIPROT ID:
O15217

Full Name:
glutathione S-transferase alpha 4

Calculated MW:
26 kDa

Purification Method:
Antigen affinity purification

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data