

For Research Use Only

AGPAT4 Polyclonal antibody, PBS Only

Catalog Number: 16562-1-PBS



Basic Information

Catalog Number:

16562-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9834

GenBank Accession Number:

BC020209

GeneID (NCBI):

56895

UNIPROT ID:

Q9NRZ5

Full Name:

1-acylglycerol-3-phosphate O-acyltransferase 4 (lysophosphatidic acid acyltransferase, delta)

Calculated MW:

378 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

AGPAT4 (Acylglycerophosphate acyltransferase 4), also known as LPAAT^δ (Lysophosphatidic acid acyltransferase delta), is a protein expressed mainly in brain and in muscle. AGPAT4 is characterized as a mitochondrial lysophosphatidic acid acyltransferase that regulates brain levels of phosphatidylcholine (PC), phosphatidylethanolamine (PE), and phosphatidylinositol (PI) (PMID: 28807933). AGPAT4 is relatively abundant in murine renal white adipose tissue (WAT) and has been linked to the regulation of TAG production, the differentiation of adipocytes, and the modulation of various pathologies including lipodystrophy, diabetes, and hepatic steatosis (PMID: 28814640).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16562-1-AP (AGPAT4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16562-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 16562-1-AP (AGPAT4 antibody) at dilution of 1:600 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 16562-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 16562-1-AP (AGPAT4 antibody) at dilution of 1:600 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 16562-1-PBS in a different storage buffer formulation.