

For Research Use Only

RPL32 Polyclonal antibody

Catalog Number: 16629-1-AP



Basic Information

Catalog Number:

16629-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9960

GenBank Accession Number:

BC070209

GeneID (NCBI):

6161

UNIPROT ID:

P62910

Full Name:

ribosomal protein L32

Calculated MW:

135 aa, 16 kDa

Observed MW:

16-18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, mouse lung tissue, HepG2 cells, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 16629-1-AP (RPL32 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.