

For Research Use Only

NEURL3 Polyclonal antibody

Catalog Number: 16648-1-AP **1 Publications**



Basic Information

Catalog Number:

16648-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10035

GenBank Accession Number:

BC012317

GeneID (NCBI):

93082

Full Name:

neuralized homolog 3 (Drosophila)

pseudogene

Calculated MW:

167 aa, 18 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

IP: mouse brain tissue,

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|---------------|-----------|-----------------------|-------------|
| Shi-Qing Zhou | 38191501 | J Exp Clin Cancer Res | WB, IHC, IF |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-NEURL3 (IP:16648-1-AP, 3ug;
Detection:16648-1-AP 1:300) with mouse brain
tissue lysate 5000ug.