

For Research Use Only

NEURL3 Polyclonal antibody, PBS Only

Catalog Number:16648-1-PBS



Basic Information

Catalog Number:

16648-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10035

GenBank Accession Number:

BC012317

GeneID (NCBI):

93082

UNIPROT ID:

Q96EH8

Full Name:

neuralized homolog 3 (Drosophila)

pseudogene

Calculated MW:

167 aa, 18 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-NEURL3 (IP:16648-1-AP, 3ug;
Detection:16648-1-AP 1:300) with mouse brain
tissue lysate 5000ug. This data was developed
using the same antibody clone with 16648-1-PBS in
a different storage buffer formulation.