

For Research Use Only

IDI2 Polyclonal antibody

Catalog Number: 16701-1-AP



Basic Information

Catalog Number:

16701-1-AP

Size:

600 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10098

GenBank Accession Number:

BC017778

GeneID (NCBI):

91734

UNIPROT ID:

Q9BXS1

Full Name:

isopentenyl-diphosphate delta
isomerase 2

Calculated MW:

227 aa, 27 kDa

Observed MW:

27-30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human skeletal muscle tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 16701-1-AP (ID12 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.