

For Research Use Only

# IDI2 Polyclonal antibody

Catalog Number: 16701-1-AP



## Basic Information

Catalog Number:

16701-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10098

GenBank Accession Number:

BC017778

GeneID (NCBI):

91734

UNIPROT ID:

Q9BXS1

Full Name:

isopentenyl-diphosphate delta  
isomerase 2

Calculated MW:

227 aa, 27 kDa

Observed MW:

27-30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human skeletal muscle tissue,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 16701-1-AP (IDI2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.