

For Research Use Only

IDI2 Polyclonal antibody, PBS Only

Catalog Number:16701-1-PBS



Basic Information

Catalog Number:

16701-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10098

GenBank Accession Number:

BC017778

GeneID (NCBI):

91734

UNIPROT ID:

Q9BXS1

Full Name:

isopentenyl-diphosphate delta
isomerase 2

Calculated MW:

227 aa, 27 kDa

Observed MW:

27-30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 16701-1-AP (ID12 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16701-1-PBS in a different storage buffer formulation.