

For Research Use Only

AFG3L1 Polyclonal antibody

Catalog Number: 16712-1-AP



Basic Information

Catalog Number: 16712-1-AP	GenBank Accession Number: BC021245	Purification Method: Antigen affinity purification
Size: 770 µg/ml	GeneID (NCBI): 172	
Source: Rabbit	Full Name: AFG3 ATPase family gene 3-like 1 (S. cerevisiae)	
Isotype: IgG	Calculated MW: 199aa, 22 kDa; 100aa, 11 kDa	
Immunogen Catalog Number: AG10134		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data