

For Research Use Only

MRFAP1L1 Polyclonal antibody

Catalog Number: 16785-1-AP



Basic Information

Catalog Number:

16785-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10217

GenBank Accession Number:

BC008087

GeneID (NCBI):

114932

UNIPROT ID:

Q96HT8

Full Name:

Morf4 family associated protein 1-like 1

Calculated MW:

127 aa, 15 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Background Information

MRFAP1L1 is thought to have a role in spermatogenesis through its interaction with TSNAX. It is involved in spermatogenesis, neuronal regulation, and genome stability.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

