

For Research Use Only

M6PR Polyclonal antibody

Catalog Number: 16795-1-AP **2 Publications**



Basic Information

Catalog Number:

16795-1-AP

Size:

70 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10267

GenBank Accession Number:

BC024206

GeneID (NCBI):

4074

UNIPROT ID:

P20645

Full Name:

mannose-6-phosphate receptor
(cation dependent)

Calculated MW:

277 aa, 31 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Anup Parchure	36173346	J Cell Biol	IF
Jian Wei	28658628	Cell Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data