

For Research Use Only

MT1E Polyclonal antibody

Catalog Number: 16831-1-AP

1 Publications



Basic Information

Catalog Number:

16831-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10440

GenBank Accession Number:

BC009699

GeneID (NCBI):

4493

UNIPROT ID:

P04732

Full Name:

metallothionein 1E

Calculated MW:

61 aa, 6 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Huan-Xin Yang	40550961	Acta Pharmacol Sin	WB,IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

