

For Research Use Only

# RP9 Polyclonal antibody, PBS Only

Catalog Number:16832-1-PBS



## Basic Information

Catalog Number:

16832-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10441

GenBank Accession Number:

BC025928

GeneID (NCBI):

6100

UNIPROT ID:

Q8TA86

Full Name:

retinitis pigmentosa 9 (autosomal dominant)

Calculated MW:

221 aa, 26 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 16832-1-AP (RP9 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16832-1-PBS in a different storage buffer formulation.



WB result of RP9 antibody (16832-1-AP; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-RP9 transfected HepG2 cells. This data was developed using the same antibody clone with 16832-1-PBS in a different storage buffer formulation.